



Casualty Care Rescue Randy

In association with Nasco Healthcare

The Casualty Care Rescue Randy is a full-body manikin that is uniquely designed for enhanced realistic training on the procedures that treat the 3 most preventable causes of death:

- 1: Massive Bleeding
- 2: Tension Pneumothorax
- 3: Airway Obstruction

- Compliant with TCCC/TECC scenarios
- Rugged, reliable, and hyper-realistic®
- 100% hand-made in the USA

- Full-body manikin with articulated joints
- Weight-distributed to mimic the human anatomy
- Holds 3-4 liters of blood
- Pumping system with regulator simulating human blood flow pressure (2-3) psi
- Manual foot operated blood pumping system allows 3 wounds to bleed simultaneously
- Optional Electric Deluxe blood pumping system allowing 4 wounds to bleed simultaneously
- Comes with blood powder (makes 5 gallons)
- Three (3) repairable neck skins
- Ten (10) repairable tracheas
- One (1) repair kit

*Please note, our products and methods may be covered by one or more of the following US Patents: 8,840,403; 9,336,693; 9,865,179; 9,916,774; 9,997,087; 10,325,524; 10,360,817; 10,360,818; 10,665,135; 10,726,743; 10,847,058, 10,977,962; 11,151,902 and US Patents Pending, as well as Foreign Patents: AU 2011271667; CA 2807472; CA 2894673; DE 60 2013 049 550.5; EP2589037 (DE); EP2589037 (FR); EP2589037 (UK); EP2852941 (FR); EP2852941 (UK).



CONTACT US DIRECTLY
(858) 244-0559



4705 RUFFIN RD
SAN DIEGO, CA
92123



WWW.STRATEGIC-OPERATIONS.COM



Casualty Care Rescue Randy:
 Life-like training solutions for
 real-life situations.

Airway Management

- Nasopharyngeal (Nose)
- Oropharyngeal (Mouth)
- Cricothyroidotomy (Throat)

Massive Hemorrhage

- Femoral and Brachial extremity wounds
- Inguinal and Carotoid junctional wounds
- Treatment Interventions
 Extremity and junctional tourniquets, wound packing and pressure bandage

Tension Pneumothorax

- Needle Decompression at the 4th/5th intercostal space in the anterior axillary line
- Needle Decompression at the 2nd/3rd Intercostal space in the mid-clavicular line

Penetrating Chest Wound

- Allows for occlusive dressing

This manikin is intended to be used as an instructional aide in the teaching and demonstrating of live rescue training. Actual product may vary slightly from photo. Nasco reserves the right to change product color, materials, supplies, or function as needed.



CONTACT US DIRECTLY
(858) 244-0559



4705 RUFFIN RD
 SAN DIEGO, CA
 92123



WWW.STRATEGIC-OPERATIONS.COM